

ELLIOTT ELECTRIC SUPPLY

We Deliver...Lower Cost, Quality Products, & Personal Service

2310 N. Stallings Dr. 75964-0000, TX Nacogdoches Phone: 936-569-7941 Fax: 936-560-4685



Manufacturers Part Number

UPC

LR200M 2" LR FM5 Conduit Body

Catalog Number	LR200M
Manufacturer	Crouse-Hinds
Description	Eaton Crouse-Hinds Series Condulet Form 5 Conduit Outlet Body, Malleable Iron, LR Shape, 2"
Weight per unit	5.3 (lbs/each)
Product Category	Bodies & Covers
eatures	
dimensions	10.5000 IN X 3.3100 IN X 4.1200 IN
Form	Form 5
laterial, Color, and Finish	
Finish	Electrogalvanized With Aluminum Acrylic Painted
imensions and Weight	
Height	3.31
Hub Size	
Length	10.50
Width	4.12
escriptions	
Description	2" LR FM5 CONDUIT BODY
extra long description	CRS-H LR200M 2 FORM 5 LR COND BODY
Features	Form 5 malleable iron conduit bodies, covers and gaskets from Eaton's Crouse-Hinds Division are used in conduit systems to act as pull outlets for conductors being installed, provide openings for making splices and taps in conductors, make 90 degree bends in conduit runs, and provide access to conductors for maintenance and future system changes. Form 5 conduit bodies are manufactured in trade sizes 1/2" to 4", and are interchangeable with Appleton Form 35 conduit bodies. They are also available with a h
Long Description	Eaton Crouse-Hinds series Condulet Form 5 conduit outlet body, Malleable iron, LR shape, 2"
Product Type	2 Form 5 LR Cond Body Malleable
Nanufacturer Information	
Brand	Eaton
GTIN	00782274756080

LR200M

782274756080



standard

ELLIOTT ELECTRIC SUPPLY

We Deliver...Lower Cost, Quality Products, & Personal Service

2310 N. Stallings Dr. 75964-0000, TX Nacogdoches Phone: 936-569-7941 Fax: 936-560-4685

Taxonomies, Classifications, and Categories	
Category Description	FORM 5 CONDULETS, COVERS, & GASKETS
Туре	Conduit outlet body
Packaging	
Carton	1
Weight Per each	5.3
Uses, Certifications, and Standards	
Application	Commercial / Institutional Buildings / Structures - Commercial /
	Institutional Buildings / Structures - Other

UL 514B, CSA C22.2, CUL